

Navina - Fecal Incontinence Insert

Description

The Fecal Incontinence Insert is a soft, flexible, liquid-filled insert designed to prevent fecal incontinence. The device is non-sterile and for single use only.

Indications for use

The Fecal Incontinence Insert is indicated for the management of Accidental Bowel Leakage (ABL) due to bowel incontinence. The rectal insert is designed for self-insertion to seal and help prevent the involuntary leakage of stool from the rectum.

Risks

- Possible mild irritation or discomfort of the anus or anal canal
- Possible migration or displacement of the insert into the anal canal or rectum
- Unexpected expulsion

Warnings, Precautions, and Contraindications

Fecal Incontinence may be medically treated. The Fecal Incontinence Insert is not intended to replace surgical options to treat accidental bowel leakage. Consult your physician for professional advice.

- The Fecal Incontinence Insert is contraindicated for use in the pediatric population.
- This device is non-sterile and for single use only. Reuse of the device could result in device breakage.
- Do not use if the package is open or damaged.
- Do not use if you have a silicone allergy.
- Do not use with suppositories.

Discontinue use of the Fecal Incontinence Insert and consult your physician if you have:

- Bloody stool or bleeding from the rectum.
- Bleeding hemorrhoids.
- Anal/rectal pain, irritation, spasms, tearing, or fissure.
- Displacement of the Fecal Incontinence Insert into the anal canal or rectum which has not been expelled with the next bowel movement.
- Undergone rectal surgery within the last 4 weeks.

Navina - Fecal Incontinence Insert

Reference	Type	Specification
6930400	Navina - Fecal Incontinence Insert, Standard	Box of 30 pieces
6930500	Navina - Fecal Incontinence Insert, Large	Box of 30 pieces
6930600	Navina - Fecal Incontinence Insert, Starter Pack	Box of 5 Standard and 5 Large pieces



At Wellspect we develop innovative continence care solutions that change people's lives. We are committed to inspire our users to build self-confidence and independence as well as good health and well-being. We have been leading the industry for over 30 years with our product brands LoFric® and Navina™. We create reliable and user-friendly products for bladder and bowel management with as little environmental impact as possible. We passionately strive to become climate neutral and work closely together with users and health-care professionals who constantly inspire us to improve our products and services in a sustainable way, now and for the future.

Wellspect. A real difference.


For more information about our products and our initiative Advancing Continence Care Together (ACCT), please visit [Wellspect.com](https://www.wellspect.com).

Join the conversation on Facebook and Instagram.

[wellspect.com](https://www.wellspect.com)



Importer: Wellspect HealthCare. A division of Dentsply IH AB, Aminogatan 1, P.O. Box 14SE-431 21 Mölndal, Sweden, Phone +46 31 376 40 00



EC/REP: MPS Medical Product Service GmbH, Borngasse 20, Braunfels 35619, Germany



Manufacturer: Minnesota Medical Technologies, 2446 Henry Road Northwest, Stewartville, MN 55976, 507-533-0366 www.mnmedtech.us